



TEMP-CAL

YSI Temperature Probe Simulator



Temp-Cal pocket size temperature simulator designed to test and validates YSI-700 compatible units. It simulates 5, 15, 25, and 40 Centigrade accurately.

It is accurate and easy to use. Connect the Temp-Cal to the instrument where the YSI 700 Probe is connected using the cable. Select the temperature on the Temp Cal. Read the YSI compatible thermometer. The reading should be +/- 1%.

110 Toledo Street, Farmingdale, NY 11735

1-800-547-6557 (US) Phone: 631-531-0100, Fax: 631-531-0101

www.NetechCorporation.Com